

CJ-210S+CJC-450 Vial Production Line



Liquid vial production processors design compliances with cGMP and Pic/s of liquid injection product guidelines. The processors are applicable for low viscosity and volume below 300mL.

The filling volume and inner plug placing parameters are computerized with access level authorization to stabilize products quality.

In order to minimize particle of pneumatic accessories polluting environment, inner plug placing procedure applies vacuum carrying and motion controlling. The capping procedure is isolated from the filling with plugging zone.

All stainless steel parts contact with product raw material are made by SUS316L and all sterilized parts are capable to be 121°C autoclave sterilized.

Specification

Model	CJ-210S Automatic Filling & Plugging machine	CJC-450 Automatic Capping Machine
Filling Volume	1ml-20ml, 20ml-50ml, 50ml-100ml,	--
	100ml-150ml, 150ml-200ml,	
	200ml-300ml	
Capacity	30-60 bottles/min	30-60 bottles/min
	(Based on filling volume and liquid characteristic)	
Power Consumption	220V single phase, 50/60HZ 2.5KW	220V single phase, 50/60HZ 1.3KW
Vacuum pressure (supplied by user)	75 NL/min, -0.6 bars	--
Compressed air (supplied by user)	120 NL/min, 4-6 bars	--
Machine Dimensions (LxWxH)	2800x1230x2193 (mm)	2200x1136x2192 (mm)
Machine Weight	950KGS	500KGS