

TCS-I Tablet Dedusting Unit



Functions

- A specialized machine is for removing the powder dust of tablets.
- Widely applicable to medicine, chemical, food industry and etc.

Features

- TCS-1 deducted is the supporting equipment for high speed tableting machine for medicine.
- With functions of high output efficiency, sealing up and low noise.
- Can be used along with air vacuum device.
- Can remove powder attached on the tablets.
- Its height is suitable for different types of tableting machines.
- Easy operation and cleaning.
- Meet PIC/S GMP standards.

Specification

MODEL	TCS-I
Max. output efficiency (tables/hr)	500,000 pcs/hr
Tablet dedusting distance	300 mm
Adjustable height of base	420 ~ 600 mm
The plane distance between feeding inlet and discharge outlet	250 mm
Machine max. diameter	320 mm
Machine max. height	950 mm
The diameter size of supporting base	315 mm
Output voltage	220/50Hz 110V/60Hz, 0.5Amp
Machine size (L x W x H)	350 x 350 x 900 m/m
Packing Size (L x W x H)	410 x 410 x 1,000 m/m
Machine total weight	30 KG